

Lentiviral Integration Sites Analysis

Advantages

- Multiplex samples for cost-effective results
- Effective workflow and fast turnaround time
- Multiple approaches to meet different goals
- Qualitative and quantitative analysis
- Comprehensive bioinformatics analysis

The Introduction of Lentiviral Integration Sites Analysis


Lentiviral vectors are the most common gene delivery vehicles used for stable genetic modification in cells. These vectors are integrated into the host genome to express therapeutic transgenes. Although this integration appears random, literature has suggested that there are in fact hotspots for lentiviral integration. To avoid deleterious effects on cellular functions (e.g., integration near oncogenes or within key genes), it is important to identify the exact integration sites of lentiviral vectors after cells have been infected. Integration site analysis is the key tool to assess the biosafety of vectors for gene therapy and the clonal tracking fate of genetically modified cells *in vivo*.

CD Genomics has developed a comprehensive set of tools for viral integration site analysis by combining PCR and sequencing-based approaches with proprietary bioinformatics pipelines. We are able to identify the position of integration sites and determine the integration frequency of viral vectors with unmatched precision and sensitivity. Through our end-to-end services, we are committed to providing the highest quality data and the most reliable analysis.


Method	Sample Requirement	Sequencing Strategy	Data Output	Application
LM-PCR	1x10 ⁶ cells or >500 ng gDNA	Illumina PE150	1 Gb	For CAR-T cell genome
Whole genome sequencing			≥150 Gb	For monoclonal cell line

Workflow

A. Schematic diagram of LM-PCR library construction




B. Bioinformatics analysis flow chart




Demo Results

Base Quality and Error Rate Distribution Plot(coh003_Data_L1)



A. Sequencing quality value and error rate distribution diagram




C. Original data composition


Index	Chr	Start	End	SupportReads	Func.refGene	Gene.refGene	GeneModel.refGene	ncbiRefSeq	cosmic79	BreakPoint	Distance2TSS	Distance2CDS
1	chr1	885774	886074	LTR4	Intergenic	NOC2L	-	-	-	LTR7339 chr1:886045	8605	283
2	chr1	915736	916046	LTR4	exonic	PERM1	-	-	-	LTR7339 chr1:915811	1379	2042
3	chr1	909090	909230	LTR4	Intergenic	ISG15L	dbp-1001 dbp-4280	-	-	LTR7339 chr1:909196	2074	1070
4	chr1	871676	871732	LTR36	Intergenic	ASPM	-	-	-	LTR7339 chr1:871717	16174	1439
5	chr1	874076	874376	LTR36	Intergenic	ASPM	-	-	-	LTR7339 chr1:874277	18868	1127
6	chr1	876821	876921	LTR4	exonic	ASPM	-	None-GHG 351	-	LTR7339 chr1:876880	21919	0
7	chr1	1048194	1048495	LTR32	Intergenic	Ctcf10B	-	-	-	LTR7339 chr1:1048391	3237	1980
8	chr1	1074231	1074631	LTR4	ncRNA_intronic	LINC01342	-	-	-	LTR7339 chr1:1074414	1825	1375
9	chr1	1086888	1086988	LTR30	Intergenic	LINC01342	dbp-1001 dbp-4280	-	-	LTR7339 chr1:1086832	14272	6244
10	chr1	11011030	11011300	LTR4	Intergenic	LINC01342	dbp-1001 dbp-4280	-	-	LTR7339 chr1:11011135	1153	447

E. Integration site detection and annotation results

coh003_Data_L1 GC content



B. Base content distribution diagram



D. Distribution plot of mapped reads in chromosomes